




integrity

MOBILITY SOLUTIONS LTD

2025

PRODUCT GUIDE





Contents

Integrity have been designing and brand manufacturing for many of the healthcare industries leading suppliers **worldwide for almost 25 years** including a number of the **UKs largest providers** to hospitals, community equipment stores and mobility dealers.

Our range of products are available to own brand by full mixed containers supplied directly from our wholly owned **manufacturing base in China.**

Call: **0800 043 0050**

sales@integritymobilitysolutions.com

www.integritymobilitysolutions.com

Around the Home



UEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90583	ADL-7118	Over Chair Table Laminated Table Top, Steel Frame, 4 Lockable Castor wheels, Tilt, Small side table can be attached on Left or Right Hand side. Black Powder Coated & Oak Laminate		900	720 - 1220	700 - 910					15	12.5		
90290	ADL-7016	Overbed Table Height Adjustable, NO Castors, Cream Powder Coated Steel Frame c/w Beech Veneer Top	1pc/Carton	697	413	630 - 930					5	7		1314
90291	ADL-7032	Overbed Table Height Adjustable, c/w 4 x 2" Braked Castors. Cream Powder Coated Steel c/w Beach Veneer Frame	1pc/Carton	697	413	680 - 980				50	5	7		1314
90072	ADL-7006	2 Tier KD Dinner Trolley	1pc/Carton	501	445	845 - 945				100	140	7		720
90454	ADL-7001	"Vinyl Footstool Height Adjustable, c/w 4 x Ferrule"	1pcs/Carton	500	315	380-530	380-535				100	4.5		910
90455	ADL-7039	"Vinyl Footstool Height Adjustable, c/w 4 x Castors"	1pcs/Carton	500	345	375-525	375-525			40	100	4.5		910





UEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90167	ADL-7003	"Perching Stool HF Padded Seat, Adjustable Height"	5pcs/Carton	415 - 460	561 - 620	561 - 705	540-686	500-646			160	3.8		950
90168	ADL-7011	"Perching Stool with HF Seat Steel Arms and No Back, Adjustable Height"	Stacked /Stacked on Pallets 11pcs/ Stack 44pcs/ Pallets	415 - 460	560 - 620	723-868	540 - 685	500-645	510		160	4.4	1144	1188
90169	ADL-7007	Perching Stool with HF Seat, Steel Arms and Steel Back - Adjustable Height	Stacked / Stacked on Pallets 10pcs/ Stack 40pcs/ Pallets	415 - 460	560 - 620	815 - 960	540 - 685	500-645	510		160	4.8	1144	1188
90117	ADL-7051	Perching Stool with HF Seat, Steel Arms and HF Back - Adjustable Height	Stacked / Stacked on Pallets 10pcs/ Stack 40pcs/ Pallets	415 - 460	560 - 620	830 - 975	540 - 685	500-645	510		160	4.7	1040	1188
90471	ADL-7119	Bariatric PU Perching Stool Steel Arms and Steel Back. Adjustable Height (290kg Weight Limit)	Stacked / Stacked on Pallets	500 - 550	690-750	950 - 1095	570 - 715	535-680	596		254	7.6	270	324
90327	ADL-7036	Bariatric Perching Stool with HF Seat, HF Back and Steel Arms. Adjustable Height (290kg Weight Limit)	Stacked / Stacked on Pallets	500 - 550	690-750	950 - 1095	570 - 715	535-680	596		254	7.6	270	324
90270	ADL-7018	Perching Stool, PU Seat, Only	5PC/Carton	415 - 460	560 - 620	545 - 685	535 - 680	500-645			160	3.5		950



UEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90271	ADL-7013	Perching Stool with PU Seat, Steel Arms and No Back - Adjustable Height	Stacked / Stacked on Pallets	415 - 460	560 - 620	723-868	535 - 680	500-645	510		160	4	1144	1188
90272	ADL-7024	"Perching Stool, PU Padded Seat, Steel Back & Arms, Adjustable Height PU Colour: RAL 7040 (Window Grey) Frame: Steel, Powder Coated Frame Colour: RAL 9003 (Signal White)"	Stacked / Stacked on Pallets	415 - 460	560 - 620	815 - 960	535 - 680	500-645	510		160	4.6	1144	1188
90273	ADL-7073	"Perching Stool, PU Padded Seat, PU Padded Back & Steel Arms, Adjustable Height PU Colour: RAL 7040 (Window Grey) Frame: Steel, Powder Coated Frame Colour: RAL 9003 (Signal White)"	Stacked / Stacked on Pallets	415 - 460	560 - 620	855 - 1000	535 - 680	500-645	510		160	4.9	1040	1188
90275	ADL-7071	"Perching Stool with PU Seat, PU Back and PU Padded Arms - Adjustable Height"	Stacked / Stacked on Pallets	415 - 460	560 - 620	855 - 1000	535 - 680	500-645	510		160	4.9	1040	1188
90346	ADL-7028	Bariatric Perching Stool with PU Seat, PU Back and Steel Arms - Adjustable Height (290kg Weight Limit)	Stacked / Stacked on Pallets	465-515	650-710	790 - 935	485-630	440-585	596		254	7.3	270	324
90141	ADL-7017	Half Step	Carton	430	735	112					190			2145



UJEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90141	ADL-7068	2 Tier Trolley with Brakes	1pc/Carton	530	595	770 - 910				150	125			695
90533	ADL-7014	Freedom 2 Bed	2pc/Wooden Case	2230	980	335 - 735		295 to frame / 145 to actuator			SWL 213 / MAX USER 178			96
90558	ADL-7080	Freedom 2 Bed Lever	1pc/Carton	126	286	713					160	3.1		2280

Bathroom Aids



90292	ADL-7022	Premier Bed Rail Adjustable Height	1PC/Carton	833	418	808					160	10		720
90354	ADL-7010	Double Bedstick Adjustable Width Powder Coated White	4pcs/Carton	335	960 - 1555	440					140	3.6		1100



UEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90355	ADL-7019	Double Bedstick Adjustable Width Chrome Plated	4pcs/Carton	335	960 - 1555	440					140	3.6		1100
90321	ADL-7041	Bedstick Slatted Bed	5PC/Carton	822	622	500					190	2.4		2375
90322	ADL-7054	Bedstick Divan Bed	5PC/Carton	805	622	646					190	2.7		2375
90324	ADL-7031	Straps (PAIR) for UEL90321/90322	Carton	2500	25						50	0.25		150000
90521	ADL-7082	Bed Grab Rail	1pc/Carton	653	435	510					165	13.1		352
90325	ADL-7049	Bed Lever c/w Straps	1pc/Carton	610	362	561					160	3.1		1665

UEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
	90520	ADL-7030	Deluxe Lightweight Aluminium Commode Adj Height c/w PU Back, PU Horseshoe Seat with Plug, PU Arms and Potty c/w Lid.	3 per Carton	495	520 - 600	850 - 1000	470 - 620	505		180	8		720
	90096	ADL-7078	Stacking Commode Chair, Height Adjustable, Fixed Arms, HF Welded Padded Seat & Back. HF Padding Colour/Material: Brown & Smooth Powdered Coated Steel. Frame Colour: Brown Tube Clip, Anti Rattle Ring & DIP Moulded Ferrule c/w M10 Steel Washer inside Colour: Black RAL 9005	Stacked / Stacked on Pallets 11pcs/Stack 44pcs/Pallet	525-545	551	785-930	452-602	430-580	507	160	6.4	924	996
	90323	ADL-7023	Stacking Commode Chair, Height Adjustable, REMOVABLE Arms, HF Welded Padded Seat & Back. Supplied c/w Commode Pan & Lid, Commode Pan to have Metal Handle HF Padding Colour/Material: Brown & Smooth Powdered Coated Steel Frame Colour: Brown Tube Clip, Anti Rattle Ring & DIP Moulded Ferrule c/w M10 Steel Washer inside Colour: Black RAL 9005	Stacked / Stacked on Pallets 11pcs/Stack 44pcs/Pallets	525 - 545	551	785 - 930	452-602	430-580	507	160	5.7	672	747
	90522	ADL-7057	Bed Bumpers (Blue with black mesh)	5 Pairs/ Carton	2000	710	10					2		1665
	90560	ADL-7034	Bed Bumpers (Cream with white mesh)	5 Pairs/ Carton	2000	710	10					2		1665
	90397	ADL-7035	Postural Wedge	6pc/ Carton	530	275	210					0.62		2496

Bathing / Showering / Toileting Aids



UEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90040	ADL-7107	6" Bath Seat	8 pcs/Carton	290	450	150	150							
90041	ADL-7038	8" Bath Seat	6 pcs/Carton	280	450	205	205				190	2.5		2106
90042	ADL-7044	12" Bath Seat	6 pcs/Carton	280	450	305	305				190	2.8		1368
90210	ADL-7106	26" Bath Board	5 per carton	240 (9.5")	660 (26")	20 (3/4")								
90211	ADL-7108	27" Bath Board	5 per carton	240 (9.5")	690 (27")	20 (3/4")								
90212	ADL-7110	28" Bath Board	5 per carton	240 (9.5")	710 (28")	20 (3/4")								
90254	ADL-7105	26" Shower Board	5 per carton	360 (14")	660 (26")	20 (3/4")								
90255	ADL-7109	27" Shower Board	5 per carton	360 (14")	660 (27")	20 (3/4")								
90256	ADL-7104	28" Shower Board	5 per carton	360 (14")	660 (28")	20 (3/4")								



UEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90115	ADL-7025	Shower Stool, Height Adjustable. Powder Coated Steel c/w PP Seat Colour: White	4pcs/Carton	370 - 395	529-540	560 - 710	480 - 630	450-600	485		127	3.3		1280
90392	ADL-7058	KD Teak bath stool. Adjustable set height: 13" - 18" Max user weight: 300 kg	White carton / colour printed label	555-610	555-610	380-510					136	2.44		2702
90348	ADL-7008	Integrity AP Aluminium Shower Chair c/w PU Back, PU Horseshoe Seat c/w Plug, 4 x Braked Nylon Castors, Fold Up Arms, Removable Back, Removable Swing Away Leg Rests c/w Nylon Footplates, Potty. Seat width: 46cm. Max user weight: 159kg	1pc/Carton	860	569	985-1010	510-535	470-495	469	125	160	13.4		451
90349	ADL-7077	Integrity SP Aluminium Shower Chair c/w PU Back, PU Horseshoe Seat c/w Plug, 2 x Braked Nylon Castors, 2 X 24" MAG Wheels, c/w PU Tyres and Q.R Pins, Fold Up Arms, Removable Back, Removable Swing Away Leg Rests c/w Nylon Footplates, Potty. Seat width: 46cm Max user weight: 159kg	1pc/Carton	1044	690	985-1010	510-535	470-495	469	125	160	13.4		216
90350	ADL-7042	Integrity AP HD Aluminium Shower Chair c/w PU Back, PU Horseshoe Seat c/w Plug, 4 x Braked Nylon Castors, Fold Up Arms, Removable Back, Removable Swing Away Leg Rests c/w Nylon Footplates, Potty. Seat width: 56cm. Max user weight: 260kg	1pc/Carton	860	669	985-1010	510-535	470-495	569	125	260	17		342
90513	ADL-7004	Integrity AP XXL Aluminium Shower Chair c/w PU Back, PU Horseshoe Seat c/w Plug, 4 x Braked Nylon Castors, Fold Up Arms, Removable Back, Removable Swing Away Leg Rests c/w Nylon Footplates, Potty. Seat width: 60cm Max user weight: 300kg	1pc/Carton	860	715	985-1010	510-535	470-495	615	125	300	19		324



UEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90351	ADL-7040	Integrity Self propel wheel kit consisting of 2 x 24" Wheels c/w PU tyres, Aluminium brakes, wheel mountings and Wheel cnversion kits	1 set / Carton	625	195	630								816
90341	ADL-7027	Toilet Frame Height Adjustable, Complete with Injection Moulded Seat	Stacked on Pallets	443-471	607-612	622-771	428 - 578	276-426	523		190	3	840	1056
90532	ADL- 7087	Toilet Frame Height Adjustable & Width Adjustable, Complete with Injection Moulded Seat	Stacked on Pallets	443-471	570-720	622-771	428 - 578	276-426	523		190	3	840	1056
90071	ADL-7012	Aluminium Swivel Bather. Adjustable Width	1pc/Carton		770 - 820	530		150	600-700		150	7.3		440
90120	ADL-7037	Toilet Frame, Clip on Seat & Chute, Height Adjustable.	Stacked / Stacked on Pallets 15pcs/ Stack	465 - 495	558	690 - 835	483-628	470-615	498		190	4		1067
90438	ADL-7015	"Toilet Surround - Adjustable Width, Adjustable Height"	6pcs/Carton	450 - 455	584 - 734	645 - 795			486-636		190	4		1080



UEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90565	ADL-7101	Toilet Surround - Adjustable Width, Adjustable Height With Floor Fixing	6pcs/Carton	450 - 455	584 - 734	645 - 795			486-636		180	4		1080
90028	ADL-7075	"Toilet Surround - Free Standing, Adjustable Height"	6pcs/Carton	450 - 455	524	645 - 795			426		180	4		1044
90564	ADL-7102	Toilet Surround - Free Standing, Adjustable Height With Floor Fixing	6pcs/Carton	450 - 455	524	645 - 795			426		180	4		1044
90116	ADL-7053	Bath Step, 4 Part	1pc/Carton	375	480	110					250	2.51		3278
90092	ADL-7047	Static Shower Chair Height Adjustable	Stacked / Stacked on Pallets	480-505	620 - 645	775-925	440-590	420-570	500		160	4	924	990
90099	ADL-7060	"Mild Steel Chrome Plated Glideabout. 4 Castors, Removable Backrest, Grey Aperture Seat with Padded Vanity Cover (Pad), Removable Arm Rests, Removable Leg Rests and Potty, Grey Vinyl Upholstery"	1pc/Carton	840	525	985	525	460	425	125	165	14		420



UEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90294	ADL-7081	Chrome Plated Height Adjustable Glideabout.	1pc/Carton	840	525	985 - 1085	525 - 625	460	425	125	160	15		523
90391	ADL-7021	"Mild Steel Grey Fleck Powder Coated Glideabout. 4 Castors, Removable Backrest, Removable Armrests, Removable Leg Rests and Potty. Black Vinyl Upholstery"	1pc/Carton	880	525	985	525	465	425	125	165	14		420
90517	ADL-7005	"Mild Steel Grey powder coated Glideabout. 4 Castors, Removable Backrest, Aperture Seat with Padded Vanity Cover (Pad), Removable Arm Rests, Removable Leg Rests and Potty, Vinyl Upholstery"	1pc/Carton	922	530	962	528	460	430	125	165	13.1		420
90595	ADL-7085	"Mild Steel Grey powder coated Glideabout. 4 Castors, Removable Backrest, Aperture Seat with Padded Vanity Cover (Pad), Removable Arm Rests, Removable Leg Rests and Potty, Vinyl Upholstery Height Adjustable"	1pc/Carton	922	530	985-1085	528	460	430	125	165	13.1		420
90430	ADL-7084	Mobile Shower Commode	1pc/Carton	519	533	780-880	450/520	400-500	451	100	160	7.5		336
90146	ADL-7063	Drop Down Rail	1pc/Carton	692	100	750					190	4		5496



UEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90519	ADL-7062	KD Glideabout	1pc/Carton	902	520	908	530	465	428	200	190	10.3		612
90398	ADL-7050	Swivel Bath Stool	5pc/Carton	380	380	410 - 580	410 - 580				135	3.30		3170
90416	ADL-7064	Aluminium Knock Down Bath Bench	4pc/Carton	370	545	516 - 641	420 - 545				160			1280
90597	ADL-7100	Male Urinal Bottle c/w lid		126	98	225								

Mobility Aids



UEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90386	ADL-7066	3 Wheeled Aluminium Rollator c/w. 200 X 32mm PU Castor and Bag (Foldable)	1pc/Carton	670	685	800 - 980			532	200	125	5.2		558
90387	ADL-7069	4 Wheeled Aluminium Rollator. c/w 200mm PU Castor, Bag, Seat and Back Rest (Foldable)	1pc/Carton	670	630	800 - 980	550	300	445	200	145	7.1		513
90518	ADL-7079	3 wheeled Mild Steel Rollator	1pc/Carton	670	715	780-930			510	200	125	6.6		558
90527	ADL-7046	"Lightweight Aluminium 4 Wheel Rollator Material: Aluminium Product weight: 6.9 kg Material of wheels: PP+PVC Finish: Metallic Size(mm): 750x650x875-975 Seat Height(mm): 545 Seat Size(mm): 460x250 Front wheel size(mm): 8" Rear wheel size(mm): 8" Max user weight (kgs): 136 KGS"	1pc/Carton	750	650	875 - 975	545			200	136	6.90		1015
90526	ADL-7074	"HD Aluminium 4 Wheel Rollator Material: Aluminium Product size: L700 x W780 x H910-1010 (mm) Seat Height: 580mm Seat Width: 560mm Shopping bag, cane holder attached Front wheel 8 inch; rear wheel 10 inch c/w TPR tyres Height adjustable handles Loop brake system for running and parking Box Dimensions: L725 x W295 x H935(mm) Net Weight: 9.6Kgs Gross weight: 12.1Kgs Max user weight: 200KGS"	1pc/Carton	700	780	910 - 1010	580			200 front 250 Rear	200	9.60		1270





UEL Part Number	ADL Part Number	Description		Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90559	ADL-7083	4 Arm Walker		1pc/Carton	795	645	1100-1225	535	530	445	200	136	9		1015
90531	ADL-7072	Wheelchair Back Pack		20pc/Carton	930	540	110						0.5		14400
UEL Part Number	ADL Part Number	Description	Colour	Method of Packing	Open Height	Folded Height	Folded Width	Folded Depth	Handle	Foot	Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised		
90586	ADL-7113	360° Walking Cane	Blue		765mm – 940mm in 25mm adjustments	380mm	130mm	87.5mm	(L) 130mm x (W) 31,5mm x (H) 25mm – 38mm	87.5mm diameter					
90587	ADL-7114		Red												
90588	ADL-7115		Black												
90589	ADL-7117		Violet												

Moving & Handling



UEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90528	ADL-7043	RS1 Patient Turner	1pc/Carton	472	475	914				75	150	11.5		728



UEL Part Number	ADL Part Number	Description	Method of Packing	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	Seat Height (mm)	UnderSeat Clearance (mm)	Width Between Arms (mm)	Castor Size (mm)	Max Weight (KG)	Net Weight (KG)	QTY in 40HQ - Palletised	QTY in 40HQ - Unpalletised
90338	ADL-7056	RS3 Patient Mover c/w Support Seat	1pc/Carton	829	700	940 - 1010	515 - 585	515 - 585		125	150	23.5		225
90380	ADL-7059	RS2 Patient Transfer Assist	1pc/Carton	635	607	1186				125 (F) / 25 (R)	160	14.5		510
90330	ADL-7076	Ladder Belt (Sliding Support Strap) for use on RS2 Patient Transfer Assist	1pc/Carton	2000	100	20						0.4		10000
90329	ADL-7033	Patient Handling Belt for use on RS2 Patient Transfer Assist 3 sizes available S,M,L	1pc/Carton	"S:1050 M:1300 L:1500"	"S:975 M:1150 L:1350"	100						"S:0.47 M:0.5 L:0.56"		"S:11250 M:9135 L:8100"
90345	ADL-7048	RS4 Patient transfer assist	1pc/Carton	820	530 closed / 730 open	1060	620 Under / 670 Over			75	185			352
90344	ADL-7086	RS5 Patient Transfer Assist	1pc/Carton	815	615	970	620 Bottom / 670 Top		423	75	185	21		330

Continuously developing new products

In an effort to supply you with the latest in innovative ADL products, our team of designers and engineers are continuously developing and expanding our product range.

Ask us about our upcoming products.

Call: 0800 043 0050

Email: sales@integritymobilitysolutions.com

Visit: www.integritymobilitysolutions.com



integrity

MOBILITY SOLUTIONS LTD

www.integritymobilitysolutions.com